Tefal

BODYSIGNAL

Tefal

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
Your appliance

The drawings are designed merely to illustrate the features of BodySignal and BodySignal Glass and do not strictly reflect reality.

The TEFAL SAS company reserves the right to modify the characteristics or components of the product at any time, in the interest of the consumer..

Description

1. Selection: Male/Female.
2. Age indicator.
3. Height indicator, in cm.
4. Display of the identification numbers for 4 users.
5. Display of the different values: weight and fat mass in kg.
6. Graphic representation of the fat mass value and recommended minimum and maximum limits.
7. Selection and scrolling keys.
8. Validation key for entering personal characteristics.
9. Key to clear the last personal characteristics entry.
10. VISIO CONTROL: indicator light for the evolution of the user's body make-up.
11. Electrodes for positioning the feet.
12. Battery charge indicator.

A new approach to weight

From early childhood, your weight is a decisive element in your physical condition, athletic performances and feeling of bodily wellness. But you cannot effectively control your weight over time if you do not know your body make-up (water, muscle, bone and fat mass).

Muscle and fat mass, which are the main components of weight, are indispensables to life in a balanced distribution which varies with sex, age and height.

- Muscles enable you to perform physical activity, lift objects and move easily. The greater your muscle mass, the more energy your muscles consume, which can help you to control your weight.
- Your fat mass works as insulation, protection and energy reserve.
It is inert and stores excess energy if you eat more than you burn.

Physical exercise blocks it, inactivity increases it. In large quantities, it becomes dangerous to your health (cardiovascular risks, cholesterol, etc.).

Sustainable weight loss

In the context of a diet, BodySignal helps you to monitor your loss of fat mass while making sure that your muscle mass is maintained. Indeed, a poorly suited diet can lead to muscle loss. But muscles are what consume energy.

If you lose muscle, you burn up fewer calories and your body then stores up fat. You therefore cannot lose fat mass sustainably without maintaining your muscle mass. BodySignal lets you take control of your diet and manage your body's balance over time.

Monitoring growth and preventing obesity in children

BodySignal is an effective way of screening for weight problems, starting in early childhood.

In fact, this period of growth is when the human body's parameters change the most, and it is often difficult to detect problems... And yet a propensity toward obesity, for example, is easier to treat if it is detected early on... keeping the child from becoming overweight in adulthood.

BodySignal is designed to monitor body make-up starting at the age of 6.

Measurement principle

When you stand barefoot on the scales, a very weak electric current is transmitted inside your body by the electrodes.

This current, running from one foot to the other, runs into difficulties when it has to cross through fat. This resistance by the body, called bioelectric Impedance, varies with the sex, age and height of the user and makes it possible to determine fat mass and limits.

There is **no danger** with this appliance, as it operates with a weak current (four 1.5 V batteries).

There are risks of erroneous measurements, however, in the following cases:

- People with fever.
- Children under the age of 6 (the appliance is not calibrated for people below this age).
- Pregnant women (although there is no danger).
- People undergoing dialysis, people suffering from oedemas in the legs, people suffering from deformity.

Caution: this appliance must not be used by people with pacemakers or other electronic medical devices due to the risks of disturbances caused by the electric current.

BODYSIGNAL's functions


BodySignal provides personalised monitoring:

- It measures and displays:

- Your weight in kg
- The fat mass in kg contained in your body.


This result is displayed in graph form for your minimum and maximum limits (calculated for your sex, age and height).

The objective is to make sure that your result falls within the limits.

You have too much fat mass.

Modify your eating habits and get regular exercise.

Your fat mass is below the lower limit.

You might want to increase your fat mass, as it constitutes an energy reserve for your body.

- It indicates the true evolution of your weight in fat and muscle mass with the "VISIO CONTROL" indicator light.

Each time you weigh yourself, the appliance memorises your body make-up and gains a better understanding of the natural variations related to your everyday life (food, hydration, rhythm of life, etc.). These "insignificant" weight variations do not change your body make-up in any way.

BodySignal eliminates these variations and warns you only if there is a real change in your fat or muscle mass (when on a diet, for example).

Thus, your weight or your body make-up can vary slightly without your needing to be warned.

- These scales store your personal characteristics (sex, age, height) in the memory and **recognise you automatically** the next time you weigh yourself..
- It can **identify 4 different users**, including **children age 6 and above**.

Interpretation of VISIO CONTROL

Losing or gaining weight can be considered as a good thing or a bad thing:

In fact, you can lose weight by losing fat mass, which is a positive development for your weight. But you can also lose weight by losing muscle, which is a negative development for your weight. You can also gain weight by gaining fat mass, which is a negative development, but you can also gain weight by gaining muscle, which is a positive development.

If your weight does not vary or if your variations are insignificant, you are considered as being stable.

• Weight loss:

either  a loss of fat mass. ➡ Your diet and/or exercise are having a healthy effect on your body make-up.

VISIO CONTROL is green

or  a loss of muscle. ➡ Modify your diet, with emphasis on compound carbohydrates and physical exercise.

VISIO CONTROL is pink

• Weight gain:

either  an increase in muscle. ➡ Your physical activity is having a healthy effect on your body make-up.

VISIO CONTROL is green

or  an increase in fat mass. ➡ Modify your diet, which is probably out of balance, emphasising fruits and vegetables, reduce your sugar and fat consumption and why not take up a sport?

VISIO CONTROL is pink

• Stability:

 stable. ➡ Your body make-up is stable. If you want to modify your body make-up, see the previous indications.

VISIO CONTROL is orange

NB : VISIO CONTROL lights up completely when the batteries are changed and after you weigh yourself for the first time to indicate that your personal data have been recorded in the memory.

Recommendations for use

To obtain complete satisfaction from this appliance, it must be **used as follows**:

- **You should always weigh yourself under the same conditions, preferably once a week, in the morning** (so as not to take into account the body's fluctuations during the day related to food, digestion or physical activity), **15 minutes after you get up**, so that the water in your body has time to be distributed in your lower limbs.
 - Do not weigh yourself after physical effort. **Wait 6 to 8 hours** for a significant weight measurement or weigh yourself before the effort.
 - **You MUST be barefoot when weighing yourself** (direct contact with the skin). Your feet must be dry and positioned properly on the electrodes.
 - **It is important to make sure that there is no contact:**
 - between your legs (thighs, knees, calves), between your feet.
- If necessary, slide a sheet of paper between your legs.

Changing your personal characteristics (age or height of children, for example)

BodySignal is switched off.*

- 1- Press once on ENTER.
- 2- Press on < or > to select your memory.
- 3- Press on ENTER for 5 seconds until the little person starts to flash.
- 4- The little person flashes: If you want to change the sex, press on < or >, then once on ENTER. Otherwise, press once on ENTER.
- 5- The age flashes: If you want to change the age, press on < or >, then once on ENTER. Otherwise, press once on ENTER.
- 6- The height flashes: If you want to change the height, press on < or >, then once on ENTER. Otherwise, press once on ENTER.

The results of your previous weighing will be redisplayed automatically, and then BodySignal switches off.


Clearing a memory

When the 0 flashes to indicate the memory number, this means that there is no more free memory space. Clear one of the 4 memories and re-enter your parameters.

BodySignal is switched off.*

- 1- Press once on ENTER.
- 2- Press on < or > to select the memory.
- 3- Now press on the CLEAR key for 5 seconds.
- 4- Clr is displayed.
- 5- All the data in the selected memory are erased.
- 6- Wait for the appliance to switch off before weighing yourself again.

** BodySignal Glass needs to be activated to access this function. For this, step on the scales, wait for the complete display of your personal data before stepping off the scales and pressing on CLEAR, then continue the process from step 2.*

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